

Typical Restriction Digest of Genomic DNA

For a total reaction volume of 40 μ l¹

1. 10 μ g genomic DNA
2. Enzyme - use 20-30 total units (keep final glycerol concentration at 5% or less)²
3. 10x reaction buffer (supplemented with 10X BSA) - usually supplied by the enzyme vendor.
4. RNase - use 1 μ l of 10 μ g/ μ l stock
5. Spermidine - use 0.4 μ l of 1 M stock (10 mM final)

Make a mix sufficient for the total number of samples (plus 10%). For example, if you have 27 samples make a mix for 30³:

example

Eco RI	2 μ l	x30 =	60 μ l	(20 units/ μ l)
10 X Enzyme Buffer	4 μ l	x30 =	120 μ l	
RNase	1 μ l	x30 =	30 μ l	
Spermidine	0.4 μ l	x30 =	12 μ l	

Mix, place 7.4 μ l in each tube bearing genomic DNA plus TE (volume of 32.6 μ l)⁴

Finger flick to mix. *DO NOT VORTEX.*

Pulse spin in microcentrifuge to collect.

Incubate 37°C 6 hrs to overnight.^{5,6}

Notes:

¹ dependent upon well size of gel comb.

² enzymes typically supplied in 50% glyceol

³ to allow for pipetting wastage

⁴ use genomic tips to pipette the DNA

⁵ temperature dependent upon endonuclease

⁶ the use of an incubator oven for the digestion reactions is preferable to use of a water bath to avoid condensation on the inner lids of the reaction tubes from altering the reaction concentrations

Entered by HKS from DD's notebook 3/24/99